

How to Make Patient Liver Health “Wallet” Cards


There are several different wallet cards available in a single file. Wallet cards should be printed double sided on card-stock (110 lb paper).

We’ve found that using brightly colored paper makes the wallet cards more visually appealing and allows patients to easily distinguish them from one another.

After you have printed the wallet cards, follow the instructions below for laminating and trimming the cards.



1 Place the double-sided sheet of wallet cards in the lamination pouch.



2 Next place the lamination pouch in the protective cardboard lamination sleeve.



3 Insert the lamination sleeve into the pre-heated lamination machine, flap-side first.



4 Remove the laminated sheet from the cardboard sleeve and allow to cool.



5 Align the crosses on the wallet card sheets with the cutting guide on the paper trimmer.



6 You will receive more uniform results if you first cut the cards vertically and then horizontally.



7 The finished wallet card size is 2” x 3.5”.



Wallet cards may be displayed in acrylic cases or packaged in sets of 12. The cases may be purchased at most office supply stores*.

*The display holder can be purchased at Corporate Express. The item number is DEF70801.